



Multimodal Acoustofluidics – Combining feedback control and ultrasound to move single cells!



Alexander Edthofer
Ph.D. Student





Thierry Baasch
Assistant Professor

Acoustics, Dynamics, and Control, Department of Biomedical Engineering





LUND
UNIVERSITY

Acoustic forces move cells to the closest potential minima





Multimodal single cell manipulation of calcein-stained cancer cell (DU-145) – First of a kind!




20 frequencies between 2.1 - 2.5 MHz. The experiment took 57 min, or 846 steps. Each step contained a 2s pulse at 5 V_{pp}

Time averaged image of cancer cell manipulation (with trajectory overlaid)




What if individual cells could be moved and manipulated?

Single cell analysis



Possibility 1: Select and extract a single cell from a sample for downstream analysis.


Tissue engineering



Possibility 2: Assemble cells into specific shapes, creating a scaffold.


This requires having precise control!

How do we control cell movement?




By applying a machine learning algorithm with feedback control!



The algorithms – Epsilon Greedy



The algorithms – Quadratic Programming




We track the particle/cell for every step that it takes



- 20 μm PS particles
- Selected trajectory
- Old particle position
- New particle position
- Target position
- Search area

Single particle manipulation of 20 μm polystyrene



30 modes between 2.1 - 2.6 MHz. The experiment took 13.51 min, or 413 steps. Each step contained a 0.5s pulse at 5 V_{pp}